medi GmbH & Co. KG
 Medicusstraße 1
 95448 Bayreuth
 Germany
 T +49 921 912-0
 F +49 921 912-780
 ortho@medi.de
 www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61-2 9890 8696
 F +61-2 9890 8439
 sales@mediaustralia.com.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: +32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 597, Rue Duvernay, Verchères
 QC Canada J0L 2R0
 T: +1 450-583-3317 / +1 800 361 3153
 F: +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

medi Danmark ApS
 Vejlegardsvej 59
 2665 Vallensbaek Strand
 Denmark
 T +45-70 25 56 10
 F +45-70 25 56 20
 kundeservice@sw.dk
 www.medit danmark.dk

medi Bayreuth Espana SL
 C/Canigo 2-6 bajos
 Hospitalet de Llobregat
 08901 Barcelona
 Spain
 T +34-932 60 04 00
 F +34-932 60 23 14
 medi@mediespana.com
 www.mediespana.com

medi France
 Z.I. Charles de Gaulle
 25, rue Henri Farman
 93297 Tremblay en France Cedex
 France
 T +33-1 48 61 76 10
 F +33-1 49 63 33 05
 infos@medi-france.com
 www.medi-france.com

medi Hungary Kft.
 Bokor u. 21.
 1037 Budapest
 Hungary
 T +36 1 371-0090
 F +36 1 371-0091
 info@medi.hu
 www.medi.hu

medi

Levamed® stabili-tri

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Important notes

This medical device **MD** is made for single patient use only. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Do not wear the product over open wounds, and use it only as instructed by your doctor or orthotist.

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Levamed® stabilis-tri

Intended Purpose

Levamed stabilis-tri is an ankle brace for mobilisation in a defined position with parts that can be removed.

Indications

All indications which require stabilisation or limitation of the range of motion of the ankle with subsequent functional mobilisation of the injured or impaired structures, such as:

- Ankle sprains (grade II to III)
- Chronic instability of the upper ankle and lower ankle
- Post-operative care following an operation on the ligamentous apparatus of the ankle

Contraindications

None known at present.

Risks / Side Effects

If the aids are too tight, this can cause pressure sores or constriction of blood vessels or nerves. Consult your doctor before using this device if you have the following conditions:

- A condition affecting the application area
- Paraesthesia, circulatory disorders
- Lymph drainage disorders

Tightly fitting aids can cause skin irritation (pressure, sweat, material).

Intended patient groups

Healthcare professionals should provide care to the adults and children, applying the available information on the measurements/sizes and necessary functions/indications and in line with the information provided by the manufacturer, and acting under their own responsibility.

1. Phase – acute phase (inflammation phase)

Immobilisation of the foot is absolutely necessary in order to counteract inflammations and marked swellings on the ankle occurring immediately after the injury. Fully equipped with stabilising element and strap system, this orthosis provides maximum support for relief and thus supports the reduction of swelling and pain.

2. Phase – repair phase (proliferation phase)

The mobilisation of the foot can be gradually increased, improving the muscular stabilisation of the ankle joint and assisting further reduction of swelling in the affected structures. By removing the stabilising elements, the ankle joint is given more freedom, permitting controlled movement. The combination of the stiff tongue and webbing continues to provide necessary support. In order not to re-injure structures that have healed during the day, it is however recommended to continue using the stiff element at night as well.

3. Phase - remodulation phase

In the final therapy phase, the webbing is removed to allow further freedom of movement. The aim is to continue to get the affected structures used to loads and strengthen them. The stiff tongue and the knitted fabric continue to support the body's own joint stabilisation proprioceptively.

Putting on the orthosis

- (Fig.1) Attach the middle of the strap to the front part of the bottom of the foot. Here one should ensure a marked trapezoid shape. This way, the strap end that is marked with „“ lies on the outer side of the foot and that the point symbols are visible.
- (Fig.2) Now fasten the stabilisation element (thermoplastically formable) into the knitted fabric and make sure that it is fixed to the ground at the two hook and loop points and at the calf brace.
- (Fig.3 and Fig.4) Grip the orthosis from above at both sides and secure the stabilising element with the knitted fabric. Ensure that the orthosis is wide open and now pull it onto the foot until the heel lies properly in the shell.
- (Fig.5) When closing the orthosis, please ensure that the tongue is under the knitted fabric, as with a conventional laced shoe. During lacing, the foot should be in a 90° position to the lower leg (depending on pain).
- (Fig.6) Now close the orthosis. To do this, click the laces into the eyelets and fix at the end with a bow. Make sure that the orthosis fits tightly to the foot, but avoid any pinching of the foot caused by pulling too tightly.
- (Fig.7 and Fig.8) Now grasp the end of the strap lying to the outside of the foot, „“. Guide this diagonally over the instep, backwards and upwards, and fasten it above the ankle on the fleece.
- (Fig.9 and Fig.10) Do the same with the second strap end „“, which is then fastened next to the first strap end.
- (Fig.11) The length of the strap can be individually adapted by cutting to size. To do this open the hook and loop fastenings on the strap ends and cut the strap accordingly.

Taking it off

- First loosen the strap and the lacing.
- In order to make taking it off easier, please make sure that the orthosis is wide open.
- Hold the tongue with one hand, with the other grasp the knitted fabric and the stabilising element.
- Now push the orthosis down over the foot.

Care instructions

Hook and loop fasteners should be closed for washing. We recommend using a laundry bag. Fabric softeners, fats, oils, lotions, ointments and soap residues can affect the material, hinder the climate comfort effect and cause skin irritations.

- Wash the product by hand, preferably using a medi clean detergent, or in delicate cycle at 30°C using a mild detergent without fabric conditioners.

- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Material composition

Polyamide, Polyester, Elastane

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

Your medi team

Wishes you a speedy recovery!

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR). The traceability of this product is assured via a UDI code .

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