medi GmbH & Co. KG
 Medicusstraße 1
 D-95448 Bayreuth
 Germany
 T +49 921 912-0
 F +49 921 912-780
 ortho@medi.de
 www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61-2 9890 8696
 F +61-2 9890 8439
 sales@mediaustralia.com.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Posthoornstraat 13/1
 3582 Koersel
 Belgium
 T: + 32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 597, Rue Duvernay, Verchères
 QC Canada J0L 2R0
 T: +1 450-583-3317 / +1 800 361 3153
 F: +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

medi Danmark ApS
 Vejlegardsvej 59
 2665 Vallensbaek Strand
 Denmark
 T +45-70 25 56 10
 F +45-70 25 56 20
 kundeservice@sw.dk
 www.meditandmark.dk

medi Bayreuth Espana SL
 C/Canigo 2-6 bajos
 Hospitalet de Llobregat
 08901 Barcelona
 Spain
 T +34-932 60 04 00
 F +34-932 60 23 14
 medi@mediespana.com
 www.mediespana.com

medi France
 Z.I. Charles de Gaulle
 25, rue Henri Farman
 93297 Tremblay en France Cedex
 France
 T +33-1 48 61 76 10
 F +33-1 49 63 33 05
 infos@medi-france.com
 www.medi-france.com

medi Hungary Kft.
 Bokor u. 21.
 1037 Budapest
 Hungary
 T +36 1 371-0090
 F +36 1 371-0091
 info@medi.hu
 www.medi.hu

medi

M.4[®] AGR

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Important notes

This medical device **MD** is made for single patient use only. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Do not wear the product over open wounds, and use it only as instructed by your doctor or orthotist.

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M.4[®] AGR

Intended purpose

M.4 AGR is a brace with a frame that is intended to reinforce and stabilise the knee joint using limitation of extension and flexion.

Indications

All indications that require physiological stabilisation and reinforcement of the knee joint and/or stress relief of the joint in at least two planes, such as:

- For genu recurvatum

Contraindications

None have been noted to date.

Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

Intended patient groups

Healthcare professionals should provide care to the adults and children, applying the available information on the

measurements/sizes and necessary functions/indications and in line with the information provided by the manufacturer, and acting under their own responsibility.

Important Note

Only change or modify the extension limits on the instruction of your healthcare provider or orthopaedic technician. To avoid hyperextension, always insert extension inserts at 0°. Furthermore, ensure the correct positioning of the brace: the central parts of the hinge should be in line with the top of the kneecap; the hinge should be positioned behind the medial lateral line of the leg.

Instructions for the orthopaedic technician:

Adjusting the extension:

Extension limitation inserts (grey; extension limits):
0°, 10°, 20°, 30°.

- 1 The medial adjustment (inner) and lateral adjustment (outer) of the degrees should always be equal. Extension inserts at 0° are already integrated into the brace.
- 2 To change the inserts, loosen the cover screws, remove the cover and unscrew the insert screws. To easily remove the inserts, gently bend the knee device.
- 3 The extension inserts are marked R-E (RIGHT) or L-E (LEFT) and indicate the respective degrees. Position the insert in such a way that the letter L or R is the same as the hinge plate. Push the

new extension insert between the two plates and fasten it with the screws.

- 4 Bend and flex the hinge to check that the inserts have been correctly positioned. Lastly, reposition the hinge cover.

Replacing the rubber washers:

The rubber washer inside the hinge can be replaced if it becomes worn or in order to increase or reduce the resistance by selecting the correct density of the mechanical device (increase or decrease, respectively, the hardness of the rubber washers). Before replacement, select the rubber washer with the most suitable hardness (shore - hardness of the rubber washer: 80 shore = min, 90 shore = medium, 95 shore = max), cross-referencing the weight and recurvatum degree as shown in Table 1.

NB: The rubber washer positioned in the medial part (inner) and the rubber washer located in the lateral part (outer) should always have the same hardness and degree of wear.



Table 1

Instructions for use:

- Remove the Velcro tab located on the internal part of the hinge.
- Unscrew the six screws and then remove the plate.
- Carefully extract the rubber washer located in the anterior part of the hinge and insert the new rubber washer as shown in the photograph.
- Reposition the hinge plate making sure to correctly fitted and fasten it using the four 4x6 screws (the shorter screws) for the upper part and the two 4x8 screws for the lower part.
- Lastly, affix the Velcro oval (provided in the packaging) to the cover.

Applying The Brace

(fig. 1) The tibia support is positioned with the wider part towards the top and the concave part cushioning the inside. Remove the paper and affix the two adhesive hook and loop tabs to the outer part of the tibia support so that they correspond with the two tibial straps (#1 and #4).

(fig. 2) Open all of the hook and loop closures. Sit on a chair and position the brace around the affected leg; flex the leg so that it forms a 45° angle. Position the hinges in such a way that the centres of the hinge are positioned in line with the top part of the kneecap.

(fig. 3+4) Keep the knee in a flexed position and support the two hinges at the back of the leg. Fasten the upper rear strap (#2) and then the upper front strap (#6).

(fig. 5+6) Keep the knee in a flexed position and support the two hinges at the back of the leg. Fasten the upper rear strap (#2) and then the upper front strap (#6).

(fig. 7+8) Fasten the lower rear strap (#5); and then the upper rear strap (#3).

For knee devices worn without a tibia support, follow the numbering of the straps: 1, 2, 3, 4, 5 and 6.

To prevent the brace from slipping, gently flex the knee and check that the Velcro straps have been correctly fastened.

Care instructions

Hook and loop fasteners should be closed for washing. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the paddings and covers by hand, preferably using medi clean washing agent.
 - Do not bleach.
 - Leave to dry naturally.
 - Do not iron.
 - Do not dry clean.
- Check regularly that the joint screws are secure, tighten them if necessary.

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Material composition

Avional, velours, polyester, PU-foam, Elastomer

Liability

The manufacturer's liability will become

void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

Your medi Team

Wishes you a quick return to full fitness.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR). The traceability of this product is assured via a UDI code [UDI](#).

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