medi GmbH & Co. KG
 Medicusstraße 1
 D-95448 Bayreuth
 Germany
 T +49 921 912-0
 F +49 921 912-780
 ortho@medi.de
 www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61-2 9890 8696
 F +61-2 9890 8439
 sales@mediaustralia.com.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: + 32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.
 Room 1703, T1, CES WEST BUND CENTER,
 No. 277, Long Lan Road
 Xuhui District
 200232 Shanghai
 People's Republic of China
 T: +86-21 50582319
 F: +86-21 50582319

UK Responsible Person
 medi UK Ltd.
 Plough Lane
 Hereford HR4 OEL
 Great Britain
 T +44-1432 37 35 00
 F +44-1432 37 35 10
 enquiries@mediuk.co.uk
 www.mediuk.co.uk

M.4[®]s OA comfort
4-Punkt Hartrahmenorthese zur Varus- / Valgusentlastung und Stabilisierung · 4-point knee orthosis for varus or valgus pressure relief and stabilisation
 Gebrauchsanweisung. Instructions for use. Mode d'emploi. Instrucciones de uso. Instruções para aplicação. Istruzioni per l'uso. Gebruiksaanwijzing. Bruksanvisning. Návod k použití. Instrukcja po użyciu. Instrukcja z wykorzystaniem. Návod na použitie. دليل الاستخدام. 使用说明. Instrucțiuni de utilizare. שימוש. הוראות שימוש. Bruksanvisning.

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Important notes
 This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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M.4s OA comfort

Intended purpose

M.4s OA comfort is a brace with a frame designed for relieving load on and stabilising the knee joint.

Indications

All indications that require load relief of a compartment via the 3-point principle, physiological support and stabilisation in accordance with the 4-point principle and, where required, limitation of the range of motion of the knee joint, such as:

- Severe medial or lateral compartment osteoarthritis (class III at least) in conjunction with complex ligament instability.
- Before HTO (corrective osteotomy) for unicompartamental load relief (brace test)
- Following cartilage surgery or cartilage regeneration
- Following surgical intervention on the meniscus

Contra-indications

Bilateral osteoarthritis of the knee

Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the

area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Putting on the brace

- Open all straps and pull the numbered ends out of the slots on the belt. Tip: To stop the straps from sticking to each other, you should fasten the hook and loop part after this step.
- Sit on the edge of a chair and pull the brace over the affected leg. Bend your leg at approximately 70°.
- Position the brace so that the centre of the joint (middle of the extension wedges) sits level with the top edge of the kneecap and slightly behind the leg's lateral central line. Please make sure that the brace has been placed on the leg straight and is not twisted.
- Fasten strap 1 first and then strap 2.
- Keep your leg bent and check again that the brace is positioned correctly.
- Then, fasten straps 3 to 6 in order.
- Take a few steps in order to ensure

that the brace is sitting correctly. If required, the straps can be adjusted, in the same order, to ensure that the brace is fitted correctly.

- To take off the brace, open all of the straps and pull the brace off.

Adjusting for Varus/Valgus Load Relief

Have your technician explain the adjustment options for the M.4s OA to you so that you are clear on them. This will enable you to find a way to manipulate the brace safely so that you can benefit from optimal support from the brace as quickly as possible.

Setting Load Relief: for Professionals

Set the load relief mechanism to neutral. This is achieved when the pressure component lies flat against the side of the knee joint and the patient can only feel slight pressure. Ensure that both brace joints are still parallel to each other. Then increase the pressure by turning all of the setting screws a quarter-turn towards "+". This is a good base setting for the first wearing period. Let the patient walk around wearing the brace for a few minutes. If necessary, adjust the setting screws again.

Setting the restrictions for flexion and extension (can only be carried out by an orthopaedic technician)

Extension restriction wedges, silver (stretch restriction): 0°, 10°, 20°, 30°, 45°
Flexion restriction wedges, black (bend restriction):

0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

Can be immobilised at:

0°, 10°, 20°, 30°, 45°

The medial (internal) and lateral (external) degree settings must always be the same.

The brace is fitted with a 0° extension wedge. Other wedges are included with the brace.

Select the required wedge (extension wedges are marked with "Ext", flexion wedges are marked with "Flex", the relevant number of degrees is also stated).

Place the wedge in the brace joint and lock it with the screw.

Move the brace joints and check that the wedges are securely in place and are well-positioned.

Important Note:

Only change the extension and flexion restrictions under the orders of the doctor treating you. To avoid hyperextension, extension wedges or 0° wedges must always be used. The wedges set the exact radius of movement for the brace. If soft tissue conditions pose issues, then it may be necessary to use a larger wedge, due to soft tissue compression, in order to limit movement of the knee to the desired radius. Ensure that the brace joints are in the correct position: The middle of the brace joints should be level with the top of the kneecap. The brace joint must be behind the lateral central line of the leg.

Care instructions

The brace cannot be worn in saltwater. If the brace is worn in chlorinated water, it should be rinsed out thoroughly with clean water afterwards. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the paddings and covers by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

**Material composition**

Avional, velours, polyester and PU-foam.

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

**Your medi team**

Wishes you a speedy recovery!

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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