medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
F +49 921 912-0
F +49 921 912-57
info@medi.de
www.medi.de

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
P +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: +32 11 24 25 60
F: +32 11 24 25 64
info@medi-belgium.be
www.medi-belgium.be

medi Canada Inc / médi Canada Inc
597, Rue Duvernay, Verchères
QC Canada J0L 2R0
T: +1 450-583-3317 / +1 800 361 3153
F: +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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medi Norway AS
Vestheimvegen 35
4250 Koppervik
Norway
T +47 52 84 45 00
F +47 52 84 45 09
post@medinorway.no
www.medinorway.no

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

www.medi.biz/anziehen

Tips on putting on and wearing compression stockings

Your mediven compression stocking is a superior quality therapeutic garment. Here are a few application and care tips to help you:

www.medi.biz/donning

Symbolerläuterungen / Explanation of the symbols / Explication des symboles / Explicación de símbolos / Explicação de símbolos / Spiegazione dei simboli / Verklaring van de symbolen / Symbolforklaring / Symbolförklaring / Vysvětlení symbolů / Objašnjenje simbola / Пояснение символов / Açıklanan simgeler / Objašnjenje symboli / Επεξηγήση συμβόλων / Szimbólumok magyarázata / Пояснення символів / Razlaga simbolov / Vysvetlenie symbolov / معاني الرموز / 符号说明 / 符號說明 / Explicare simboluri / הסבר על הסמלים / Forklaring av symboler / Обяснение на символите / Sūmbolite selgitus / Symbolien selitykset / Simbolu skaidrojums / Simbolių paaiškinimas

Medizinprodukt / Medical product / Produit médical / Producto médico / Dispositivo médico / Prodotto sanitario / Medisch product / Medicinsk produkt / Medicinteknisk produkt / Zdravotní výrobek / Medicinski proizvod / Медицинское изделие / Tibbi ürün / Produkt medyczny / Ιατροτεχνολογικό προϊόν / Orvostechnikai eszköz / Медицинский продукт / Medicinski pripomoček / Zdravotnícky produkt / منتج طبي / 医疗器械 / 醫療產品 / Produks medis / מכשיר רפואי / Medicinsk produkt / Медицинско изделие / Meditsiiniseade / Lääkinnällinen tuote / Mediciniska ierice / Medicinos priemonė

Eindeutige Produktidentifizierung / Unique product identification / Identification unique du produit / Identificación clara del producto / Identificação única do produto / Identificazione unica del prodotto / Duidelijke productidentificatie / Entydig produktidentificering / Entydig produktidentificering / Jednoznačná identifikace výrobku / Jednoznačná identifikacija proizvoda / Однозначная идентификация изделия / Kesin ürün tanımlama / Jednoznaczna identyfikacja produktu / Σαφής αναγνώριση προϊόντος / Egyértelmű termékazonosítás / Чітка ідентифікація продукту / Enoznačna identifikacija izdelka / Jednoznačná identifikácia produktu / تحديد واضح للمنتج / 医疗器械唯一标识 / 明確的產品標識 / Identificare univocă a produsului / הזהרת "חיד" / Entydig produktidentifikasjon / Унікална ідентифікація на изделието / Unikaalne tootetähistus / Tuotteen yksiselitteinen tunnistus / Viennozzimiga produkta identifikacija / Aiškus priemonės identifikavimas

Mehrmalige Anwendung durch einen Patienten (s. auch 6. Trage- und Nutzungsdauer) / Multiple uses by one patient (see also 6. Length of time that compression clothing should be worn and used for) / Plusieurs utilisations par un patient (voir également 6. Durée de port et d'utilisation) / Uso repetido por parte de un paciente (véase también 6. Duración y tiempo de uso) / Reutilização por um único paciente (ver também a secção 6. Período de utilização e vida útil) / Utilizzo ripetuto da parte di un paziente (vedere anche Paragrafo 6 Durata dell'utilizzo) / Meermalig gebruik door een patiënt (zie ook 6. Draag- en gebruiksduur) / Centagen brug ved samme patient (se også 6. Bære- og anvendelsestid) / Flergångsbruk av en patient (se även 6. Användningstid) / Vícenásobné použití jedním pacientem (viz také 6. Délka nošení a používání) / Višekratna upotreba jednog pacijenta (v. i 6. Nošenje i trajanje upotrebe) / Многократное использование одним пациентом (см. также раздел 6 «Продолжительность ношения и использования») / Bir hasta tarafından çok kez kullanım (ayrıca bakın 6. Giyme ve kullanım süresi) / Wielokrotne stosowanie przez jednego pacjenta (patrz także 6. Okres noszenia i użytkowania) / Επανελημμένη χρήση από έναν ασθενή (βλ. επίσης 6. Διάρκεια χρήσης) / Egy beteg általi többszöri használat (lásd még: 6. Viselési és használati időtartam) / Багаторазове використання одним пацієнтом (див. також 6. Час носіння та використання) / Večkratna uporaba na enem pacientu (gl. tudi 6.

medi Polska Sp. z o.o.
ul. Łąbiedzka 22
44-121 Gliwice
Poland
T: +48-32 230 60 21
F: +48-32 202 87 56
info@medi-polska.pl
www.medi-polska.pl

medi Bayreuth Unipessoal Lda.
Rua do Centro Cultural, no. 43
1700-106 Lisboa
Portugal
T: +351 21 843 71 60
F: +351 21 847 08 33
medi.portugal@medibayreuth.pt
www.medi.pt

medi RUS LLC
Business Center NEO GEO
Butlerova Street 17
117342 Moscow
Russia
T: +7-495 374 04 56
F: +7-495 374 04 56
info@medirus.ru
www.medirus.ru

medi Trading (Shanghai) Co. Ltd.
Room 1703, T1,
CES WEST BUND CENTER,
No. 277, Long Lan Road
Xuhui District
200232 Shanghai
People's Republic of China
T: +86-21 50582319
F: +86-21 50582319

medi Sweden AB
Box 6034
192 06 Sollentuna
Sweden
P: +46 8 96 97 98
F: +46 8 626 68 70
info@medi.se
www.medi.se

Medi Turk Ortopedi Medikal İc ve Dis Tic. AS
Mustafa Kemal Mah. 2156 Sokak No: 16 / 4-5
06510 Çankaya Ankara
Turkey
T: +90 312 435 20 26
F: +90 312 434 17 67
info@medi-turk.com
www.medi-turk.com

medi UK Ltd.
Plough Lane
Hereford HR4 0EL
Great Britain
T: +44 1432 37 35 00
F: +44 1432 37 35 10
enquiries@mediuk.co.uk
www.mediuk.co.uk

medi Australia Pty Ltd
83 Fennell Street
North Parramatta NSW 2151
Australia
P: +61-2-9890 8696
F: +61-2-9890 8439
sales@mediaustralia.com.au
www.mediaustralia.com.au

medi Ukraine LLC
Evhena Sverstiuka str, 11.
Kiev 02002
Ukraine
T: +380 44 591 11 63
F: +380 44 392 73 73
info@medi.ua
http://medi.ua

medi USA LP
6481 Franz Warner Parkway
Whitsett, N.C. 27377-3000
USA
P: +1 336 4 49 44 40
F: +1 888 5 70 45 54
info@mediusa.com
www.mediusa.com

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mediven® Instructions for use

Dear Customer,
Thank you for placing your trust in us! We're delighted that you've chosen medical compression garments from medi. To ensure that your compression therapy is effective, please take a moment to read through these instructions for use carefully. Ensure that you store them carefully so that you can also refer back to the important information at a later point in time. Should you have any questions, please consult the doctor treating you, the medical professionals that you have sought care from or the medi Service team (see reverse).

1. Intended use - what is a medical compression sleeve or glove used for?

Round- or flat-knitted medical compression garment used for compression of the upper extremities, mainly for the treatment of disorders of the venous or lymphatic systems.

2. Features - how does medical compression garments work?

Medical compression garments works using a graduated pressure pattern (decreasing from the end of your limb towards your heart) to exert compression on your limbs.

3. Indications - which conditions is medical compression garments used to treat?

You have received compression garments from medi following your diagnosis. Please follow your doctor's orders and/or the orders of your medical specialist in order to obtain optimal therapeutic results. Please note that each item of compression garments has different indications and applications due to their differing material properties, how they were manufactured and the compression classes that they are available in.

In general, your mediven compression garments can be used for the following indications:

Chronic venous conditions

- Venous malformations

Thromboembolic venous diseases

- Superficial vein thrombosis
- Venous thrombosis of the upper extremity
- A post-thrombotic condition
- Post-thrombotic syndrome
- Thrombotic prophylaxis in mobile patients

Oedema

- Lymphoedema
- Post-traumatic oedema
- Post-operative oedema
- Post-operative reperfusion oedema
- Cyclic idiopathic oedema
- Lipoedema
- Conditions involving stasis due to immobility (paresis and partial paresis of the extremity)
- Oedema caused by taking medication which cannot be altered (no viable alternative)

Other indications

- Conditions following burn injuries
- Scar treatment

Please observe the following instructions:

Due to the different manufacturing methods used for medical compression garments and the different modes of action that they elicit, medi gene-

rally recommends that circular-knit compression (garments without seams) be used for venous conditions and flat-knit compression (garments with seams) be used for lymphatic system conditions. However, flat knit clothing may still be an appropriate therapy for a venous condition, for example (e.g. if the limb is very large at point and narrow at another or if skin folds are very deep). Note that this decision can only be made by a doctor.

Patient-dependent factors in particular, such as the patient's weight, the type and severity of the oedema and the nature of their connective tissue all play a role in this decision.

Therefore, our recommendation is as follows:

- The higher the patient's weight,
- The more they tend towards oedema,
- The more severe the disease,
- The softer their connective tissue,
→ The more compressive the garments material needs to be!

Your doctor or the medical professionals treating you will support and advise you in selecting the correct medical compression garments for you. This ensures that your individual needs are met, which increases your wellbeing and provides you with the best chances of therapeutic success.

4. Contraindications - are there instances in which I should not wear the medical compression garments?

Medical compression garments may not be worn if the patient has the following conditions:

- Advanced peripheral arterial occlusive disease (if one of these parameters is present: ABPI < 0.5, ankle arterial pressure < 60 mmHg, toe pressure < 30mmHg or TcPO₂ < 20 mmHg on dorsum of foot). When non-elastic materials are implemented, the patient can still try out compression stocking if they have ankle arterial pressure between 50 and 60 mmHg and they are under close clinical supervision.
- Decompensated heart failure (NYHA III + IV)
- Septic thrombophlebitis
- Phlegmasia cerulea dolens

In the following instances, the therapeutic decision should be made by weighing up the benefits and risks and by selecting the most suitable compressive equipment:

- Severe weeping dermatosis

- Intolerance to compressive material
- Severe paraesthesia in the limbs
- Advanced peripheral neuropathy (e.g. in diabetes mellitus)
- Primary chronic polyarthritis

If you are not sure whether one or more of these statements applies to you, speak to your doctor and/or the medical personnel treating you.

The following risks and the following side effects may occur:

Medical compression garments may incur the following, especially if handled incorrectly:

- Skin necrosis and
- Pressure damage to peripheral nerves.

If your skin is sensitive, the compressive equipment may cause itchiness, peeling or signs of inflammation. Adequate skincare under the compression clothing is vital. Therefore, please take particular care to observe our important notes and the instructions for putting the compression garments on (sections 7 and 8).

If the following symptoms are experienced then the compression garments should be removed immediately and clinical advice sought: If the toes are blue or white in colour, if you experience paraesthesia or numbness, increasing pain, shortness of breath and break out in sweats or acute limitation in terms of your movement

5. Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

6. How long should my mediven compression garments be worn for? How long can it be used for?

Unless otherwise specified by the doctor, wear your mediven compression stocking from morning to night for optimal therapeutic results. A

precisely graduated compression pattern (high pressure at bottom, decreasing towards the top) is required to ensure that compression garments functions optimally. Daily wearing and washing of the garments can cause the compression and elasticity of your garments, which are required to make them medically effective, to deteriorate over time. The stockings should therefore only be used for a maximum of six months. If the medical compression garments has been used for six months and a repeat prescription is issued, the patient's dimensions must be checked again by the medical specialist retailer as they could have changed due to the patient's condition or their individual circumstances. Observe our instructions for putting on the armsleeve and instructions for care to ensure that the compressive effect is maintained as well as possible (sections 8 and 10). Please don't hesitate to get in touch with us. We can give you information on the options for ensuring long-term effective therapy for venous and lymphatic system disorders using medical compression garments.

7. Instructions for use - what should I know?

- The compression garments should be issued by trained specialists with a briefing on how to handle it and how to care for the compression garments (see sections 8 and 10 on this).
- If pain or increased irritation of the skin occurs while you are wearing the compression garments, remove it immediately and consult your doctor or the medical specialists.
- Appropriate skincare must be maintained during compression therapy. This is essential. Creams or salves that are particularly greasy and also soap residue can cause skin irritation and material wear and can impair the efficacy of the compression garments. Therefore, medi also offers a skincare range (medi day, medi night, medi soft foam lotion) that is specially formulated for compression clothing. Please visit a medical specialist retailer to discuss which product is right for you.
- Compression garments with a silicone top band can cause skin irritation in people with sensitive skin (caused by perspiration in combination with mechanical stress). In order to prevent this, it is recommended that the top band is repositioned slightly several times throughout the day. If the patient uses unsuitable skincare (see also previous point) or has a lot of body hair, then the

compression garments could slip. Please visit a medical specialist retailer to discuss this, should it occur.

- Sharp fingernails and jewellery can cause damage to the knit fabric. Do not attempt to repair your compression garments yourself or to have it repaired by an external service provider, as this will render the guarantee void, without exception. Furthermore, this can impair the quality, safety and efficacy of the product. You can obtain information on correct repair by medi - if possible - from your medical specialist retailer or our services.
- In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).
- The sewn-in fabric label is an important component in ensuring that the identification and traceability of the product is guaranteed.
- Therefore, ensure that it is not removed from the compression garments under any circumstances. Do not tear it out or cut it out, as this renders any claim to the guarantee, repair or exchange null and void.
- The RAL quality mark symbol shown on the packaging of your standard product applies exclusively to our sleeves.
- Medical compression wear which has been tailored to the patient's own dimensions cannot be exchanged.

8. Care instructions - what do I need to do when putting on my mediven compression garments?

Explanatory videos with tips and tricks for putting on the compression garments correctly can be found here:

Ideally, you should put the compression garments on immediately after you wake up. If you have a shower in the morning and then apply skin care products (see "Section 10 - Skin care products" here), then you should wait a while until your skin is completely dry. This will make it far easier to put on the compression garments. Please also ensure that you have your medi fabric gloves or rubber

gloves ready. These are available from specialist medical retailers. You will use them later. They make it easier to distribute the fabric over the arm, or you can use them to put the sleeve on, right from the beginning. The gloves proffer excellent purchase and make it easier to put on the sleeve.

Tip: Ask about our handy donning aid, the medi Arm Butler.

When you are putting the sleeve on, ensure that your fingernails have been filed to a smooth, rounded shape so that the knit fabric does not become damaged. Only put your jewellery or watch on once you have put the compression garments on. Instructions for putting the compressive garments on using the medi Arm Butler

Instructions for PUTTING the compressive garments on using the medi Arm Butler

Tip: Place a non-slip surface under the medi Butler (e.g. a sponge cloth or a rubberised sheet).

A. Place the medi Arm Butler on a table or similar surface.

B. If the compression garments has a seam, roll the sleeve over the semicircular cylinder. Ensure that the seam is not central and pointing directly to you but is instead slightly offset.

- Compression sleeves for the right arm - slightly to the outer right (at approx. 4 o'clock)
- Compression sleeves for the left arm - slightly to the outer left (at approx. 8 o'clock)

C. Now lie the medi Arm Butler flat on the table. To stabilise it, place the front bracket on the front edge of the table.

D. Put your hand through the opening. Ensure that you are holding down the Butler firmly and carefully by its cuff with your other hand.

E. Insert your arm into the compression sleeve, keeping the angle the same and without turning it.
F. Check that the compression sleeve is seated correctly again. If you wish, fold up the top band or fastening.

9. Storage instructions - how long can I store my mediven compression garments, and how should I store it?

- The products must be stored in a dry place and protected from direct sunlight.
- The hourglass symbol shown on the box label or on the measurement chart shows the latest date that the compression clothing can be used

up until. For standard products the expiry date has been calculated by adding the maximum storage duration for standard products (36 months) and the maximum amount of time that the stocking can be used for (6 months).

- Models that have been tailored for the patient and made-to-measure versions are intended for immediate use. Therefore, the expiry date is calculated by adding the production, transportation and delivery duration (one month) and the maximum amount of time that the stocking can be used for (six months).

10. Care instructions - what do I need to take note of in order to maintain my mediven compression garments optimally?

- Wash your compression garments every day after you have worn it. We recommend that you use medi clean or a specially formulated hand wash detergent to ensure gentle care. Alternatively, you can use washing detergent for delicates to wash the garments by hand or in the machine (delicates cycle at 40°C) without optical bleach or conditioner.

Tip: A net washing bag offers extra protection for the knit fabric.

- We recommend that you turn your compression garments inside out when washing it. Wash it separately or with clothing that is similar in colour.
- If the compression garments is stained, we recommend that you use medi spot ex. Please do not dry clean mediven compression garments or clean it with household agents.
- After washing, you can roll the compression garments up in a towel and press it to dry it. Do not leave the damp compression garments in the wet towel, instead let it dry in the air. Do not place it in sunlight or on a radiator/heating appliance.
- Alternatively, you could also dry your mediven compression garments in your tumble dryer using the delicates cycle.
- Stretch compression garments with a seam (flat knitted garments) back into shape after washing.

Wash up to 40°C (delicates cycle)

Do not bleach

Tumble dry at low temperatures (delicates cycle)

 Do not iron

 Do not dry-clean

11. Material composition - which ingredients does my mediven compression garments contain?

- All mediven compression garments contains polyamide and elastane.
- You can find more detailed information on the fabric label sewn into your compression clothing.
- Latex free: This product does not contain any rubber that has been obtained from the latex tapped from tropical plants.
- mediven Harmony and mediven harmony SL contain lanolin. Lanolin is a natural wool wax used in skin care products."

12. Disposal - how should I dispose of my mediven compression garments?

- Please dispose of your medical compression sleeves via your general waste or household waste. There are no specific criteria for disposing of the stockings.

13. Explanation of the symbols

- Medical product
- Unique product identification
- Multiple uses by one patient (see also 6. Length of time that compression clothing should be worn and used for)
- Cannot be exchanged
- Store in a dry place
- Keep out of sunlight

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