
mediven[®] thrombexin[®] 18
mediven[®] thrombexin[®] 21
mediven[®] struva[®] 23
mediven[®] struva[®] 35


Gebrauchsanweisung. Instructions for use. Инструкции за употреба. Brugsvejledning. Kasutusjuhend. Käyttöohje. Mode d'emploi. Οδηγίες χρήσης. 使用説明書. Istruzioni per l'uso. Uputa za uporabu. Lietošanas instrukcija. Naudojimo instrukcija. 使用说明. Gebruiksaanwijzing. Bruksanvisning. Instrukcja zakłádania. Manual de instrucções. Instrucţiuni de utilizare. Инструкция по использованию. Bruksanvisning. Návod na použitie. Navodila za uporabo. Instrucciones de uso. Návod k používání. Kullanna kilavuzu. تعليمات الاستعمال. Інструкція з використання. Használati útmutató. הוראות שימוש.

1

MD


Medizinprodukt / Medical product / Медицинское изделие / Medicinsk produkt / Meditsiiniseade / Lääkinnällinen tuote / Produit médical / Ιατροτεχνολογικό προϊόν / 醫療產品 / Presidio medicale / Medicinski proizvod / Mediciniska ierice / Medicinos priemonė / 医疗器械 / Medisch product / Medicinsk produkt / Produkt medyczny / Dispositivo médico / Produs medical / Медицинское Изделие / Mediceinteknisk produkt / Zdravotnícky produkt / Medicinski pripomoček / Producto médico / Zdravotnícky prostriedek / Tibbi ürün / Медицинний виріб / منتج طبي / Orvosteknikai eszköz / מכשיר רפואי

2

UDI


Eindeutige Produktidentifizierung / Unique product identification / Унікальна ідентифікація на изделието / Entydig produktidentificering / Unikaalne tootetähistus / Tuotteen yksiselitteinen tunniste / Identification unique du produit / Στοιχός αναγνώρισης προϊόντος / 明確的產品標識 / Identificazione unica del prodotto / Jednoznačna identifikacija proizvoda / Viennozīmīga produkta identifikācija / Aiškus priemonės identifikavimas / 医疗器械唯一标识 / Duidelijke productidentificatie / Entydig produktidentifikasjon / Jednoznaczna identyfikacja produktu / Identificação única de dispositivo / Identificare univocă a produsului / Унікальний Ідентифікатор Изделий / Entydig produktidentifering / Jednoznačná identifikácia produktu / Enoznačna identifikacija izdelka / Identificación clara del producto / Jedinečná identifikace produktu / Kesin ürün tanımlama / Унікальний ідентифікатор продукту / تحديد واضح للمنتج / Egyértelmű termékazonosítás / מזהה מוצר ייחודי


3


Anzahl der Wiederverwendungen minus eins / Number of times product can be reused minus one / Брой на повторните употреби минус едно / Antal gange gentaget brug minus en / Korduvkasutuste arv miinus üks / Uudelleenkäyttöjen lukumäärä miinus yksi / Nombre de réutilisations moins un / Αριθμός εκ νέου χρήσεων μείον ένα / 重複使用次數減一 / Numero di riutilizzi meno uno / Broj ponovnih korištenja minus jedan / Atkārtoto izmantojumu skaits mīnus viens / Panaudojimo kartų skaičius atėmus vieną / 复用次数减1 / Aantal keren dat een product kan worden hergebruikt min één / Antall gjenbruk minus ett / Liczba określająca, ile razy można założyć pończochę, pomniejszona o jeden / Número de vezes que o produto pode ser reutilizado menos um / Număr de reutilizări minus una / Количество повторных циклов „использование-стирка” минус один / Antal gånger som produkten kan återvändas minus en / Počet opätovných použití minus jedna / Številov ponovnih uporab, zmanjšano za ena / Número de repeticiones de uso menos uno / Počet opakovaných použití minus jedna / Yeniden kullanim sayısı eksi bir / Кількість разів повторного використання продукту, мінус один / عدد مرات إعادة الاستخدام ناقص واحد / מספר השימושים החוזרים פחות אחד / Újrafelhasználások száma mínusz egy

4

5

6

7

DEUTSCH

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mediven® thrombexin® 18 und mediven® thrombexin® 21:

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Intended purpose

mediven® thrombexin® 18 and mediven® thrombexin® 21:

Round knitted clinical compression stocking for compression of the lower extremities, mainly for prevention of diseases of the venous system.

mediven® struva® 23 and mediven® struva® 35:

Round knitted clinical compression stocking for compression of the lower extremities, mainly for prevention and treatment of diseases of the venous system and general swelling.

Product features

These products work using a graduated pressure pattern (decreasing from the end of your limb towards your heart) to exert compression on your limbs.

Indications

mediven® thrombexin® 18:

- Physical thrombosis prophylaxis in immobile patients - pre-, intra- and post-operative

- Prevention and therapy of post-operative and post-traumatic oedema

mediven® thrombexin® 21:

- Physical thrombosis prophylaxis in immobile / partially mobile patients - pre-, intra- and post-operative
- Physical thrombosis prophylaxis in immobile patients with an increased risk of thrombosis (e.g. current venous disease, history of deep vein thrombosis)
- Prevention and therapy of post-operative and post-traumatic oedema and haematoma and for reduction of concomitant pain

mediven® struva® 23:

- Physical thrombosis prophylaxis in partially mobile patients - pre-, intra- and post-operative
- Prevention and therapy of post-operative and post-traumatic oedema and haematoma and for reduction of concomitant pain (e.g. orthopaedic or venous surgery)
- Therapy of venous leg ulcer in immobile / partially mobile patients

mediven® struva® 35:

- Physical thrombosis prophylaxis in mobile patients - pre-, intra- and post-operative
- Prevention and therapy of post-operative and post-traumatic oedema and haematoma and for reduction of concomitant pain (e.g. orthopaedic or venous surgery)
- Therapy of venous leg ulcer in partially mobile / mobile patients
- Compressive treatment of acute deep vein thrombosis

Contraindications

Advanced peripheral arterial occlusive disease (if one of these parameters is present: ABPI < 0.5, ankle arterial pressure < 60 mmHg, toe pressure < 30 mmHg or TcPO₂ < 20 mmHg on dorsum of foot), Decompensated heart failure (NYHA III + IV), Septic thrombophlebitis, Phlegmasia cerulea dolens, Massive leg oedema

Take especial care if the patient has

Severe weeping dermatosis, Intolerance to compressive material, Severe paraesthesia in the limbs, Advanced peripheral neuropathy (e.g. in diabetes mellitus), Primary chronic polyarthritis

Risks and side effects

The products may incur the following, especially if handled incorrectly: skin necrosis and pressure damage to peripheral nerves. If your skin is sensitive, the compressive equipment may cause itchiness, peeling or signs of inflammation. Adequate skincare under the compression stocking is therefore vital. Please comply with the instructions for use and inspect your skin each day, checking that your circulation is optimal. If you experience the following symptoms, the compression stocking should be removed immediately and clinical advice should be sought: Change in toe colour to blue or white, new paraesthesia or numbness, an increase in pain, shortness of breath, sudden sweating or acute limitation of mobility.

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that the persons supporting carer has received a briefing on this from healthcare professionals. Healthcare professionals should provide care to adults and children, applying the available information on measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Instructions for use

Optimal care can only be provided if the measurements are taken and applied correctly. **mediven® thrombexin® 18 and mediven® thrombexin® 21:**

The circumference of the ankle (cB) dictates the size that must be selected. For the knee stocking, the measurement of the calf circumference (cC) confirms the correct size. For thigh stockings, the selection of the size is confirmed by the measurement of the thigh circumference (cG), or this may indicate that an extra-wide stocking is needed. In the interests of better handling and legibility, measurements for garments are given in whole centimetres. Any potential gaps that may result from this are covered by measurement and production tolerances. The internationally recognised rules for rounding apply for the measurements (round down for numbers up to 4/10, round up from 5/10). This was taken into account when designing the paper measuring tapes provided by the manufacturer.

mediven® struva® 23 and mediven® struva® 35: The circumference of the ankle (cB) determines the size that must be selected. For knee stockings, the measurement of the circumference under the knee (cD) confirms the correct size. For thigh stockings, the selection of the size is confirmed by the under-the-knee circumference (cD) and the measurement of the thigh circumference (cG).

The products may only be used on intact skin or if an appropriate wound dressing is in place, and are suitable for use by multiple patients, provided that they are cleaned in advance. Contact with chemical substances (e.g. benzene, methanol, salicylic acids in high concentrations) may impair the compressive function of the products. Oils and ointments should not be applied to the limb when compression stockings are used. Compression stockings with a silicone top band can cause skin irritation in people with sensitive skin (caused by perspiration in combination with mechanical stress). In order to prevent this, it is recommended that the top band is repositioned slightly several times throughout the day. If you experience pain or increased skin irritation, take off the compression stocking immediately and consult your doctor or carer.

Length of time that compression clothing should be worn and used for

Unless otherwise stipulated by your doctor, you can wear the compression stocking day and night in order to elicit optimal therapeutic results. The compression stockings should be changed as often as necessary. Under normal circumstances, this should be every 3 days, at the least. The compression stockings can be laundered a maximum of 10 times.

Instructions for care

Launder according to the requirements of the relevant Medical Device Directive and the relevant national requirements for laundering medical devices by regarding the following pre-conditions:

- Wash at temperatures of up to 95°C
- Do not bleach
- Dry at temperatures of up to 80°C
- Do not iron
- Do not dry-clean

Damaged compression stockings must be removed from circulation and can no longer be used on patients for the set indications.

Storage

The products must be stored in a dry place and out of direct sunlight.

In combination with other products

A combination of mediven® struva® 23 and mediven® struva® 35 with mediven® struva® fix may be a good solution for holding wound dressing material in place, but this can only be executed with the approval of the doctor treating you.

Disposal

Please ensure proper and suitable disposal of the used products to prevent any risk to human beings and the environment.

The products do not contain any latex.

Complaints

In the event of any complaints regarding the product, such as damage to the fabric or a fault in the fit, please contact the manufacturer. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer and the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of Regulation (EU) 2017/745 (MDR).

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medi RUS LLC