**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**F +49 921 912-783**  
**ortho@medi.de**  
**www.medi.de**

medi Australia Pty Ltd  
 83 Fennell Street  
 North Parramatta NSW 2151  
 Australia  
 T +61-2 9890 8696  
 F +61-2 9890 8439  
 sales@mediaustralia.com.au  
 www.mediaustralia.com.au

medi Austria GmbH  
 Adamgasse 16/7  
 6020 Innsbruck  
 Austria  
 T +43 512 57 95 15  
 F +43 512 57 95 15 45  
 vertrieb@medi-austria.at  
 www.medi-austria.at

medi Belgium NV  
 Staatsbaan 77/0099  
 3945 Ham  
 Belgium  
 T: + 32-11 24 25 60  
 F: +32-11 24 25 64  
 info@medibelgium.be  
 www.medibelgium.be

medi Canada Inc / médi Canada Inc  
 597, Rue Duvernay, Verchères  
 QC Canada J0L 2R0  
 T: +1 450-583-3317 / +1 800 361 3153  
 F: +1 888-583-6827  
 service@medicanada.ca  
 www.medicanada.ca

MAXIS a.s.,  
 medi group company  
 Slezská 2127/13  
 120 00 Prague 2  
 Czech Republic  
 T: +420 571 633 510  
 F: +420 571 616 271  
 info@maxis-medica.com  
 www.maxis-medica.com

medi Danmark ApS  
 Vejlegardsvej 59  
 2665 Vallensbaek Strand  
 Denmark  
 T +45-70 25 56 10  
 F +45-70 25 56 20  
 kundeservice@sw.dk  
 www.medit danmark.dk

medi Bayreuth Espana SL  
 C/Canigo 2-6 bajos  
 Hospitalet de Llobregat  
 08901 Barcelona  
 Spain  
 T +34-932 60 04 00  
 F +34-932 60 23 14  
 medi@mediespana.com  
 www.mediespana.com

medi France  
 Z.I. Charles de Gaulle  
 25, rue Henri Farman  
 93297 Tremblay en France Cedex  
 France  
 T +33-1 48 61 76 10  
 F +33-1 49 63 33 05  
 infos@medi-france.com  
 www.medi-france.com

medi Hungary Kft.  
 Bokor u. 21.  
 1037 Budapest  
 Hungary  
 T +36 1371-0090  
 F +36 1371-0091  
 info@medi.hu  
 www.medi.hu

medi Japan K.K.  
 5-1 Nihonbashi Kabutocho Chuoku  
 103-0026 Tokyo  
 Japan  
 T: +81 3 6778 2590  
 F: +81 3 5847 7901  
 info@medi-japan.jp  
 www.medi-japan.co.jp

medi Nederland BV  
 Heusing 5  
 4817 ZB Breda  
 The Netherlands  
 T +31-76 57 22 555  
 F +31-76 57 22 565  
 info@medi.nl  
 www.medi.nl

medi Polska Sp. z o.o.  
 ul. Łąbiedzka 22  
 44-121 Gliwice  
 Poland  
 T: +48-32 230 60 21  
 F: +48-32 202 87 56  
 info@medi-polska.pl  
 www.medi-polska.pl

medi Bayreuth Unipessoal, Lda  
 Rua do Centro Cultural, no. 43  
 1700-106 Lisbon  
 Portugal  
 T +351-21 843 71 60  
 F +351-21 847 08 33  
 medi.portugal@medibayreuth.pt  
 www.medi.pt

medi RUS LLC  
 Business Center NEO GEO  
 Butlerova Street 17  
 117342 Moscow  
 Russia  
 T +7-495 374 04 56  
 F +7-495 374 04 56  
 info@medirus.ru  
 www.medirus.ru

medi Sweden AB  
 Box 6034  
 192 06 Sollentuna  
 Sweden  
 T +46 8 96 97 98  
 F +46 8 626 68 70  
 info@medi.se  
 www.medi.se

medi Trading (Shanghai) Co. Ltd.  
 Room 1703, T1, CES WEST BUND CENTER,  
 No. 277, Long Lan Road  
 Xuhui District  
 200232 Shanghai  
 People's Republic of China  
 T: +86-21 50582319  
 F: +86-21 50582319

Medi Türk Ortopedi Medikal İc ve Dis Tic. AS  
 Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5  
 06510 Çankaya Ankara  
 Turkey  
 T: +90 312 435 20 26  
 F: +90 312 434 17 67  
 info@medi-turk.com  
 www.medi-turk.com

medi UK Ltd.  
 Plough Lane  
 Hereford HR4 OEL  
 Great Britain  
 T +44-1432 37 35 00  
 F +44-1432 37 35 10  
 enquiries@mediuk.co.uk  
 www.mediuk.co.uk

medi Ukraine LLC  
 Evhena Sverstiuka str, 11.  
 Kiev 02002  
 Ukraine  
 T: +380 44 591 11 63  
 F: +380 44 392 73 73  
 info@medi.ua  
 http://medi.ua

medi USA L.P.  
 6481 Franz Warner Parkway  
 Whitsett, N.C. 27377-3000  
 USA  
 T +1-336 4 49 44 40  
 F +1-888 5 70 45 54  
 info@mediusa.com  
 www.mediusa.com

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### Important notes

This medical device **MD** is made for single patient use only. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Do not wear the product over open wounds, and use it only as instructed by your doctor or orthotist.

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Aluminium, Polyamid, Polyester, Elasthan, Elastomer

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# Spinomed®

## Instructions for the orthopedic technician

### Intended purpose

Spinomed is a brace designed to actively relieve load and correct the lumbar spine/thoracic spine in the sagittal plane.

### Indications

All indications requiring a relief and/or correction of the lumbar/thoracic spine as well as movement restriction in sagittal plane, such as:

- Osteoporotic collapse of thoracic and/or lumbar vertebrae
- Juvenile Scheuermann's disease
- Kyphosis with chronic back pain

### Contraindications

None known to date

### Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

### Intended patient groups

Healthcare professionals should provide care to the adults and children, applying the available information on the measurements/sizes and necessary functions/indications and in line with the information provided by the manufacturer, and acting under their own responsibility.

### Fitting the orthosis (this should be carried out only by an orthopedic technician!)

#### Wear the orthosis over your clothes.

The back brace does not need to be removed from the pocket during adjustment. Hold the back brace against the patient's back so that the lower edge of the brace is level with the coccyx. The upper edge should reach to approx. 3–5 cm below C7. The back brace is formable when cold and can be shaped without a tool. Start by shaping at the lordosis and then work upwards.

**Attention:** The back support must be adjusted so that it fits the patient exactly. There should be no gaps nor pressure points.

### Adjusting the straps

Let the patient put the brace on and close the abdominal pad. Now adjust the band straps, starting with the pelvic belt and moving on to the shoulder and then the elastic middle belt (Fig. 1). For each strap, undo the Velcro, pull the strap tight, and secure it again. As soon as all the straps have been correctly adjusted, they can be shortened using a pair of scissors.

### Pelvic strap

Ensure that the pelvic strap runs below the iliac crest. This strap should be tight.

### Shoulder strap

The shoulder strap should gently hold the patient erect, but should not cut in. Select the tension so that the shoulder strap lies flat against the shoulder. However, it should still be possible to push the flat hand between the strap and the shoulder. Stand to the side of the patient. It is important that you first tighten the complete strap, from bottom to the top. Loosen the strap in the shoulder area on one side of the body (Fig. 2). Now pass the strap through the tunnel element towards the back until the strap is taut (Fig. 3). Shorten the strap (Fig. 4). Now guide the strap to the front at a 180-degree-angle and fasten it (Fig. 5).

Check that the shoulder pads are seated correctly. They should be positioned so that the pad lies across the front of the shoulder area (Fig. 6).

The band strap should roll in underneath the armpit in order to prevent friction.

A flexible spiral rod has been incorporated into each shoulder strap to make it easier to put on.

These may be removed if necessary.

**Attention:** To prevent slipping of the back support, the pelvic belt must be done up tighter than the shoulder straps. If the support slips upwards sometimes, this has no impact on the function of the brace. However, the patient should pull the back support down again.

### Elastic middle strap

Adjust the middle strap to the optimum length so that the reverse buckle is located behind the armpit and the patient's shoulders are adequately held backwards.

### Checking the fit of the brace

- Back pad adjusted to fit the patient's spinal curvature

- The pelvic strap runs firmly below the hip bone
- The shoulder strap is flat against the body and the flat hand can still be inserted between the strap and the shoulder
- The shoulder pad is in the optimum position
- The elastic connecting strap is adjusted to the correct length

## Spinomed®

### Dear patient!

**Your doctor has recommended that you wear a Spinomed. In order to make application easier, please read the instructions carefully.**

### Putting on the orthosis

Wear the Spinomed over your clothes. To make it easier for you to put on the brace, it is helpful to stand in front of a long mirror.

### Putting on the Spinomed is like putting on a rucksack

- Undo the large hook and loop fastening at the front.
- Put on the orthosis as you would normally put on a rucksack (Fig. 7). The lighter side of the back brace goes against your back (Fig. 10).
- Slip your right hand into the loop on the right side of the fastening, grip the thumb opening on the left side of the fastening with the other hand and close the stomach fastening (Fig. 8). Alternatively, you can also use the thumb opening on the right side of the fastening. The brace is fitted correctly when the lower edge of the fastening is directly above the groin (Fig. 9).

### Now look in the mirror to check

- The fit of the back brace: this should

nestle closely against your spine. The bottom edge of the brace is level with the buttock cleft (fig. 10).

- The layout of the belt straps: ensure that the straps are not twisted, though they may be rotated when they are fed through the reversing clasp (Fig. 1).
- The fit of the shoulder pad; the strap must not cut into the underarm area (fig. 6).
- The reverse buckles; these must not press in at the side.

### Caution

To avoid the possibility of the back brace slipping, the pelvic strap must be tightened more than the shoulder strap. Slight upward slippage of the brace will not have any effect on the function of the orthosis. For the sake of comfort, however, you should pull the brace down again.

### Note on wearing

The Spinomed should only be worn for 1 to 2 hours per day in the habituation phase and in order to gradually restore the muscle activity.

The length of time that the brace is worn can be increased in stages over the course of the therapy.

The brace should be taken off during longer periods of sitting or lying down. If required, the brace must be realigned by your orthopaedic technician.

Your acceptance and cooperation play an important role in this form of therapy.

Therefore, should any issues occur, please discuss them with your doctor or therapeutic technician.

They may be able to make any potentially required changes to your Spinomed brace.

To ensure that this therapy is successful, you must follow your doctor's instructions at all times.

### Care instructions

The straps, shoulder pad and the large hook-and-loop fastening can be washed by hand in water with a mild detergent. The back support and pads can be cleaned using a moistened cloth.

Hook and loop fasteners should be closed for washing. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

### Material composition

Aluminum, polyamide, polyester, elastane, elastomer

### Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

### Disposal

The product can be disposed of in the domestic waste.

### Your medi team

wishes you a successful course of therapy!

In the event of any complaints regarding the product such as damage to the fabric

or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR). The traceability of this product is assured via a UDI code [UDI](#).

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